

Sonalleve MR-HIFU Clinical Summary Report

Patient Information

Name: Pig568218HIFULeftFemur ID: Pig568218HIFULeftFemur Date of Birth: 6/24/1985

 Date of Birth:
 6/24/1985

 Weight:
 20.0 kg

 Body Temperature:
 33.0 °C

Overall Conclusions

Therapy Statistics

 Treatment Date:
 6/24/2015

 Field Strength:
 3.0 T

 System Identification:
 60000

 Software Version:
 3.2.740.2311

 Report Time Zone:
 UTC-04:00

 From the First Sonication to the
 68 min 40 sec

Last:

Cell Type	Completed Sonications	Aborted Sonications
2 mm Regular Cell	6	0

Date: 6/25/2015

Signature:



Sonalleve MR-HIFU Clinical Summary Report

Treatment Information

Volume Name: Volume:	PTV 1 0.15	ml				
Cell	f[kHz]	Time	P[W]	t[s]	Test	Sonication Ending
2 mm Regular Cell 1	1200.1	12:59:58 PM 12:59:58 PM		20.0 Scan Protocol	X	Completed
2 mm Regular Cell 1	1200.1	1:07:55 PM 1:07:55 PM I	50 Non-default S	20.0 Scan Protocol		Completed
2 mm Regular Cell 2	1200.1			20.0 Scan Protocol		Completed
2 mm Regular Cell 3	1200.1	1:32:16 PM 1:32:16 PM I	60 Non-default S	20.0 Scan Protocol		Completed
2 mm Regular Cell 4	1200.1	1:56:52 PM 1:56:52 PM I		20.0 Scan Protocol		Completed
2 mm Regular Cell 5	1200.1	2:08:18 PM 2:08:18 PM I	50 Non-default S	20.0 can Protocol		Completed

Abbreviations:

Power Duration Frequenc




Sonalleve MR-HIFU Clinical Summary Report

Treatment Order

Time	Volume Name	Cell	Test	Sonication Ending
12:59:58 PM	PTV 1	2 mm Regular Cell 1	Χ	Completed
1:07:55 PM	PTV 1	2 mm Regular Cell 1		Completed
1:20:30 PM	PTV 1	2 mm Regular Cell 2		Completed
1:32:16 PM	PTV 1	2 mm Regular Cell 3		Completed
1:56:52 PM	PTV 1	2 mm Regular Cell 4		Completed
2:08:18 PM	PTV 1	2 mm Regular Cell 5		Completed




Images




Date: 6/24/2015 Time: 12:39:58 PM

T1W_FFE Son1_50W_20s_Tmax89p5




Date: 6/24/2015 Time: 1:08:18 PM

TherapySonication 2 Sonication1_coronal_50W_20s_85p7



Date: 6/24/2015 Time: 1:08:18 PM PTV 1


TherapySonication 2 Sonication1_axial_50W_75p2



Date: 6/24/2015 Time: 1:20:51 PM PTV 1

TherapySonication 3 Sonication2_coronal_80W_76p4






Date: 6/24/2015 Time: 1:20:51 PM

PTV 1


TherapySonication 3 Sonication2_axial_80W_83p0



Date: 6/24/2015 Time: 1:20:51 PM

PTV 1


TherapySonication 3 Sonication2_80W_sagittal_90p7



Date: 6/24/2015 Time: 1:32:38 PM

PTV 1

TherapySonication 4 Sonication3_60W_axial_72p2




Date: 6/24/2015 Time: 1:32:38 PM

PTV 1

TherapySonication 4 Sonication3_coronal_60W_73p4






Date: 6/24/2015 Time: 1:32:38 PM

PTV 1


TherapySonication 4 Sonication3_sagittal_76p5



Date: 6/24/2015 Time: 1:57:17 PM

PTV 1


TherapySonication 5 Sonication4_50W_coronal_61p5



Date: 6/24/2015 Time: 1:57:17 PM

PTV 1

TherapySonication 5 Sonication4_axial_50W_92p4




Date: 6/24/2015 Time: 1:57:17 PM

PTV 1

TherapySonication 5 Sonication4_50W_sagittal_72p4




Sonalleve MR-HIFU Clinical Summary Report



Date: 6/24/2015 Time: 2:08:42 PM

PTV 1


TherapySonication 6 Sonication5_coronal_50W_69p0



Date: 6/24/2015 Time: 2:08:42 PM

PTV 1

TherapySonication 6 Sonication4_axial_50W_83p7



Date: 6/24/2015 Time: 2:08:42 PM


PTV 1

TherapySonication 6 Sonication5_sagittal_50W_57p5




Date: 6/24/2015 Time: 12:39:58 PM

T1W_FFE Sonication1_coronal_4cm



Date: 6/24/2015 Time: 12:39:58 PM

T1W_FFE Sonication1_axial_4mm



Date: 6/24/2015 Time: 12:39:58 PM

T1W_FFE Sonication1_sagittal_SpotWidth_1cm



Date: 6/24/2015 Time: 12:39:58 PM


T1W_FFE Sonication2_2.5cm



Date: 6/24/2015 Time: 12:39:58 PM


T1W_FFE Sonication3_coronal_6p5cm






Date: 6/24/2015 Time: 12:39:58 PM

T1W_FFE sonication3_sagittal_1cm



Date: 6/24/2015 Time: 12:39:58 PM

T1W_FFE Sonication4_axial_SpotWidth_7mm



Date: 6/24/2015 Time: 12:39:58 PM

T1W_FFE

Sonication5_axial_SpotWidth_6mm